

NARK® II CONTROLLED SUBSTANCE ANALYSIS

Presumptive Drug Test Kits



SIRCHIE®

FINGER PRINT LABORATORIES


NARK® II Progressive System of Drug Identification

NARK® II has the capability of presumptively identifying several families of substances suspected of being abused Drugs. Designed to function as a transportable narcotics laboratory, it is available for use wherever the need for its capability might arise. Each of the tests is comprised of one or more chemical reagents. When a predictable color or series of colors occur within a specific testing sequence, a positive confirmation may be presumed.

Presumptive identification is generally recognized within our legal system as a component of probable cause. There is no drug identification system presently in use which completely eliminates the occurrence of false positives and false negatives. A forensic laboratory is required to qualitatively identify an unknown substance.

NARK® II Progressive System of Drug Identification allows testing of the following substances:

- Amphetamines
- Barbiturates
- Benactyzine
- Cannabis (Marijuana)
- Cocaine
- Codeine
- DMT
- Ephedrine
- Fentanyl
- GHB
- Heroin
- Ketamine
- LSD
- MDA
- MDMA (Ecstasy)
- Mescaline
- Methadone
- Methamphetamine
- Methaqualone
- Methcathinone
- Methedrine
- MMDA
- Morphine
- Oxycodone
- PCP
- Psilocybin (Magic Mushrooms)
- Rohypnol
- Talwin
- Valium



No. NARK200MK Kit and Components

| CAT. NO. | KIT DESCRIPTION | PRICE |
|------------|--|----------|
| NARK200MK | Master Kit, 130 Tests + Neutralizer: 2 boxes ea. test 01, 05; 1 box ea. 02, 03, 04, 07, 08, 10, 11, 12, 15 | \$192.00 |
| NARK200IK | International Kit, 130 Tests + Neutralizer: 1 box ea. test 01, 02, 03, 04, 05, 07, 08, 10, 11, 12, 13, 14, 15 | \$198.00 |
| NARK200CD | Customs/Douanes Kit, 50 Tests + Neutralizer: 1 box ea. test 01, 04, 07, 10, 15 | \$88.00 |
| NARK200MEG | Meg Kit, 50 Tests + Neutralizer: 1 box ea. test 01, 04, 07, 11, 15 | \$90.00 |
| NARK200PK | Port-A-Kit, 16 Tests + Neutralizer: 2 pouches ea. test 01, 02, 03, 04, 05, 07, 11, 15 | \$39.00 |
| CAT. NO. | ACCESSORY DESCRIPTION | PRICE |
| NARK20016 | Progressive Identification Card | \$2.75 |
| NARK20017 | Progressive Identification Wall Chart | \$6.75 |
| NARK20018 | Drug Field Test Results Log, 25-sheet pad | \$6.00 |

NARK® II CONTROLLED SUBSTANCE ANALYSIS

Popular NARK® II Field Test Reagents

NARK® II COCAINE ID SWIPES:

Presumptively identify cocaine from trace amounts of residue present on surfaces that have been in contact with cocaine. Results are immediate. Swipes have an unlimited shelf life ensuring reliable results.

Remove swipe from pouch.



Wipe suspected surface.

Presence of cocaine presumed if swipe turns blue.



| | | |
|---------|---------------------------|---------|
| NS20150 | Cocaine ID Swipes, 50 ea. | \$35.00 |
|---------|---------------------------|---------|



NARK2001
Marquis Reagent
General Screening



NARK2005
Duquenois-Levine Reagent
Marijuana



NARK2007
Scott Reagent Modified
Cocaine HCl/Crack



NARK20015
Methamphetamine Reagent
MDMA

| CAT. NO. | TEST | DESCRIPTION | PRESUMPTIVE TEST FOR | PRICE |
|-----------|------|--|--------------------------------------|---------|
| NARK2001 | 01 | Marquis Reagent, 10/box | Opium Alkaloids | \$10.40 |
| NARK2002 | 02 | Nitric Acid Reagent, 10/box | Differentiate betw. Heroin, Morphine | \$10.40 |
| NARK2003 | 03 | Dille-Koppanyi Reagent, 10/box | Barbiturates | \$13.85 |
| NARK2004 | 04 | Ehrlich's Reagent, 10/box | LSD | \$16.80 |
| NARK2005 | 05 | Duquenois-Levine, 10/box | Marijuana, Hashish, Hash-Oil, THC | \$17.70 |
| NARK2006 | 06 | Acid Neutralizer, 1 oz. (39ml) bottle | | \$2.40 |
| NARK2007 | 07 | Scott Reagent (modified), 10/box | Cocaine HCl, Cocaine Base (crack) | \$17.70 |
| NARK2008 | 08 | Methadone, 10/box | Methadone | \$15.00 |
| NARK2009 | 09 | PCP, Methaqualone, 10/box | PCP & Methaqualone | \$15.00 |
| NARK20010 | 10 | Special Opiates, 10/box | Heroin, Fentanyl, Oxycodone | \$15.40 |
| NARK20011 | 11 | Mecke's Modified Reagent, 10/box | All Heroin | \$15.90 |
| NARK20012 | 12 | Frohdes Reagent, 10/box | Talwin, Pentazocine | \$13.95 |
| NARK20013 | 13 | Chen's Reagent, 10/box | Ephedrine, Pseudoephedrine | \$16.75 |
| NARK20014 | 14 | Valium, Rohypnol, Ketamine, 10/box | Valium, Rohypnol, Ketamine | \$16.75 |
| NARK20015 | 15 | Sodium Nitroprusside Reagent, 10/box | Methamphetamine & MDMA | \$17.70 |
| NARK20019 | 19 | Mayer's Reagent, 10/box | Narcotic Alkaloids | \$10.40 |
| NARK20020 | 20 | KN Reagent (Fast Blue B Salts), 10/box | Marijuana (green plant & seeds) | \$16.75 |
| NARK20021 | 21 | GHB Reagent, 5/box | GHB | \$17.70 |