

NO. CNA1000

PORTABLE CYANOACRYLATE FUMING CHAMBER



FAST & EASY SET-UP
for On-The-Spot
Fuming at the Scene!

The CNA1000 Portable Cyanoacrylate Fuming Chamber provides over 1.2 cubic feet of fuming space to accommodate a large number of evidential items. The clear acrylic interlocking panels are laser engraved to guarantee correct orientation and are secured together with unique locking keys to ensure a tight and fast assembly. The professional design of this portable unit permits the use of any SIRCHIE cyanoacrylate fuming products (i.e., THE FINDER™ packets, OMEGA-PRINT™ compounds, CYANO-SHOT™ and the CYANOWAND™). A special access port on the front panel accommodates the CYANOWAND™ (sold separately) and is easily accessed by removing its magnetic cover.



Locking Key

CYANO-BLOC™ treatment pads are provided for maximum protection against the build-up of cyanoacrylate on the interior surfaces. With proper application, little maintenance is needed besides the occasional wipe-down of interior surfaces and reapplication of CYANO-BLOC™.



No. CNA1000 Components

CNA1000 CONTENTS:

- 6- Interlocking, Engraved, Precision Laser-Cut Acrylic Panels
- 12- Locking Keys
- 2- Hanging Rods
- 6- FR1006 Evidence Clips
- 2- CNA110A CYANO-BLOC™ Pre- and Post-Treatment Pads
- 1- CNA2000 THE FINDER™ Cyanoacrylate Pks., 5 ea.
- 1- CNA1000C Lightweight Aluminum & Black ABS Carrying Case w/snap closure;
- Dimensions: 20" x 14" x 2.125" (50.8cm x 35.6cm x 5.4cm); Weight: 9.4 lbs. (4.3kg)

CAT. NO.	DESCRIPTION	PRICE
CNA1000	Portable Cyanoacrylate Fuming Chamber	\$179.50

Copyright© 2008 by SIRCHIE Finger Print Laboratories, Youngsville, N.C. All Rights Reserved.

FLR6.08-599ENG-REV1

Call **800-356-7311** Today or Order
On-Line at **www.sirchie.com**

SIRCHIE®
FINGER PRINT LABORATORIES

FUMING PRODUCTS

CYANOACRYLATE



CNA103 & CNA102



CNA104



CNA106



CNA3006



CNA4000



CNA2000



CNA110



SCW100 Kit and Components

SCW100 KIT CONTENTS:

- 1-SCW101 CYANOWAND™ Heat Generator w/Cover
- 1-SCW103 Butane Fuel Can, 1.48 oz. (42g)
- 1-SCW104 Cartridge Extraction Tool
- 1-SCW200 CYANOWAND™ Cartridges (10 pk.)
- 1-SCW105 Carrying Case, Copolymer w/Foam Inserts; Dimensions: 10.3725" x 4.625" x 1.75" (26.3cm x 11.7cm x 4.4cm) Weight: 1 lb. (453.6g)

OMEGA-PRINT™ Products

CNA102 FUMING COMPOUND (20g) is economically ideal for departments that don't need large quantities of cyanoacrylate. Its resealable drop-dispensing cap ensures solution integrity while in storage.

CNA103 FINGERPRINT DEVELOPER for departments that process a large number of latent prints with cyanoacrylate. This 16 oz. (454g) bottle is supplied with four applicator bottles to maintain the integrity of the bulk container and to permit drop application to dispersal pads and disposable fuming trays.

CNA104 DISPERSAL PADS are for use with any cyanoacrylate liquid fuming compounds. They are not treated with sodium hydroxide and do not create highly toxic fumes. Pads packaged in quantities of 100.

CNA106 DISPOSABLE FUMING TRAYS are for use with or without dispersal pads. These aluminum trays are designed for fuming with cyanoacrylate liquid fuming compounds and are packaged in quantities of 25.

CYANO-SHOT™ Products

CNA3006 LATENT PRINT DEVELOPMENT SYSTEM utilizes cyanoacrylate to create an abundant amount of heat, which decreases fuming time while providing quality ridge detail.

CNA4000 ULTRA CYANO-SHOT™ features a fuming time lasting up to 20 minutes allowing you to develop latent prints in larger enclosed areas.

THE FINDER™ Packets

CNA2000 is a unique cyanoacrylate gel formulation sandwiched between two layers of aluminum foil. Simply peel apart to activate the fumes. Packets can be used repeatedly until exhausted (about 15 hours).

CYANO-BLOC™ Treatment Pads

CNA110 is specially formulated to protect new, unused fuming chambers (polycarbonate material or glass) from cyanoacrylate induced film.

CYANOWAND™ Kit

SCW100 features a self-igniting, butane-powered heat tool. It utilizes disposable cyanoacrylate cartridges for fuming latent prints. Simply attach a cartridge, ignite and begin fuming.

CAT. NO.	DESCRIPTION	PRICE
CNA102	OMEGA-PRINT™ Fuming Compound, 20g	\$4.96
CNA103	OMEGA-PRINT™ Fingerprint Developer, 16 oz. (454g)	\$49.50
CNA104	OMEGA-PRINT™ Dispersal Pads, 100 ea.	\$10.95
CNA106	OMEGA-PRINT™ Disposable Fuming Trays, 25 ea.	\$8.25
CNA3006	CYANO-SHOT™ Latent Print Development System, 6 pk.	\$21.95
CNA3012	CYANO-SHOT™ Latent Print Development System, 12 pk.	\$40.95
CNA4000	ULTRA CYANO-SHOT™	\$11.25
CNA2000	THE FINDER™ Cyanoacrylate Packets, 5 ea.	\$15.96
CNA2001	THE FINDER™ Cyanoacrylate Packets, 20 ea.	\$55.80
CNA110	CYANO-BLOC™ Treatment Pads, 5 ea.	\$25.00
SCW100	CYANOWAND™ Kit	\$153.95